

# LogBook

2024\_07\_22

## Sample: #1 Cu-6Sm

**Loading:** 2.7 mg of sample in the capillary

**Sample:** 60%wt CuO, 40%wt 6Sm (calcinated at 500°C, 4%at. Zr, 6%at Sm, 90%at Ce)

**Temperature calibration:** To reach a sample temp. “T<sub>sam</sub>”, you need to set a temperature “T<sub>sp</sub>”

$$T_{sp} = (T_{sam} - 7.27) / 0.772$$

For T<sub>sam</sub>= 200 °C; set T<sub>sp</sub>= 249.7

For T<sub>sam</sub>= 220 °C; set T<sub>sp</sub>= 275.7

For T<sub>sam</sub>= 240 °C; set T<sub>sp</sub>= 301.6

For T<sub>sam</sub>= 260 °C; set T<sub>sp</sub>= 327.5

For T<sub>sam</sub>= 280 °C; set T<sub>sp</sub>= 353.4

For T<sub>sam</sub>= 300 °C; set T<sub>sp</sub>= 379.3

For T<sub>sam</sub>= 320 °C; set T<sub>sp</sub>= 405.2

For T<sub>sam</sub>= 340 °C; set T<sub>sp</sub>= 431.2

For T<sub>sam</sub>= 360 °C; set T<sub>sp</sub>= 457.1

For T<sub>sam</sub>= 380 °C; set T<sub>sp</sub>= 483.0

For T<sub>sam</sub>= 400 °C; set T<sub>sp</sub>= 508.9

For T<sub>sam</sub>= 550 °C; set T<sub>sp</sub>= 703.3

For T<sub>sam</sub>= 600 °C; set T<sub>sp</sub>= 768.1

**ROI1:** Ce L<sub>3</sub>

**ROI2:** Sm L<sub>3</sub>

**ROI3:** Cu K<sub>α</sub> (this is irrelevant)

**ROI4:** Zr K<sub>α</sub>

**ROI5:**  $Fe K_{\alpha}$

**Plot PFY (cps):** ROIX/mcaLt (Lifetime of MCA detector)

**Plot Transmission:** Diode1/AS\_IC (average of the Ionch) for Cu K-edge mostly

**Plot Ref:** ionch1

### Probable errors:

- When you stop scan with CTRL+C, afterwards do “monomot3”, “mv monoh 25” and go to the energy again “mv monoE edge”

### Export these files

## 1 – Pristine RT (“001\_Cu-6Sm\_Pristine\_RT”)

**Temperature:** 28°C

**Flow:** set-point at 5 mL/min Ar (after 3am, 0 mL/min)

**Pressure:** atmospheric pressure

#	Type	Scan info	Temp( $T_{sp}$ )	Pres (bar)	Comments
1	Eiger	18keV	RT	1	20 sec
2	Ref	Mo K-edge	RT	1	Mo reference (filter3 pos 14)
3-5	EXAFS	Zr K-edge	RT	1	Mode 4, ROI 4 (BEAMSTOP IN)
6-9	EXAFS	Zr K-edge	RT	1	Mode 4, ROI 4 (mz=-1.5, LOW SIGNAL)
10-12	EXAFS	Zr K-edge	RT	1	Mode 4, ROI 4 (mz=0, NO SIGNAL)
14	Eiger	18keV	RT	1	mz=0 (NO SIGNAL, not in the sample?)
15	Eiger	18keV	RT	1	mz=-1.5,10 sec (in the sample)
16	Ref	Cu K-edge	RT	1	
17	Ref	Cu K-edge	RT	1	filter3 16
18-21	NEXAFS	Cu K-edge	RT	1	Low intensity
22-23	NEXAFS	Cu K-edge	RT	1	Low intensity (to see the quality)
24	alignment	cr2ropi	RT	1	At Cu K-edge
28	Ref	Cu K-edge	RT	1	With new values of cr2ropi
29-30	Ref	Cu K-edge	RT	1	With new values of cr2roll
End of shift (3:00 am)					
Start of shift (7:39 am, 24/07/2024))					
32	Eiger	18keV	RT	1	10sec
33	Ref	Cu K-edge	RT	1	
34-35	NEXAFS	Cu K-edge	RT	1	Fluor. probably saturated (No filt.)
36	NEXAFS	Cu K-edge	RT	1	Filter Cu 0.034 mm
40	Ref	Fe K-edge	RT	1	
41-42	NEXAFS	Ce $L_3$ -edge	RT	1	
43-44	NEXAFS	Ce $L_3$ -edge	RT	1	After changing fluor. detector

45	Ref	Cu K-edge	RT	1	
46-47	NEXAFS	Cu K-edge	RT	1	Good transmission, fluor. signal is saturating
48-49	NEXAFS	Cu K-edge	RT	1	1 sec acq. (Fluor. signal is saturating)
50-51	NEXAFS	Cu K-edge	RT	1	Al 0.06 mm (Fluor. signal is saturating)
52-53	NEXAFS	Cu K-edge	RT	1	Mo 0.01 mm (Fluor. signal is saturating)
54-55	NEXAFS	Cu K-edge	RT	1	Al 0.2 mm Filter3=9 (Fluor. signal is good)
<b>Cu K-edge measurement: Transmission is good w/o filters, Fluor. is good with Al 0.2 mm filter</b>					
77	Ref	Fe K-edge	RT	1	After realignment (Ivo magic)
78	NEXAFS	Ce L <sub>3</sub> -edge	RT	1	Mode 2, 2 seconds per point / ABORT
80-85	NEXAFS	Ce L <sub>3</sub> -edge	RT	1	Mode 2, 2 seconds per point
86-87	alignment	Sm L <sub>3</sub> -edge	RT	1	
88-93	NEXAFS	Sm L <sub>3</sub> -edge	RT	1	Mode 2, 2 seconds per point
94	Ref	Cu K-edge	RT	1	Abosorbes too much, not useful
95	Ref	Cu K-edge	RT	1	Filter3 = 22 (Cu 0.01 mm) NEW!
97-98	Ref Cu(0)	Cu K-edge	RT	1	Filter3 = 22 (Cu 0.01 mm), good reference
99-110	EXAFS	Cu K-edge	RT	1	Use Mode 4, Transmission (diode1)
111	Ref	Mo K-edge	RT	1	Energy calibration at 20keV
112-117	XANES	Zr K-edge	RT	1	Mode 2, 2 sec/point. Trans + Fluor.
118-120	alignment		RT	1	
122	Eiger	18keV	RT	1	10 sec
124	Eiger	18keV	RT	1	2 sec
125	Eiger	18keV	RT	1	60 sec
128	Ref	Cu K-edge	RT	1	Filter3 = 22 (Cu 0.01 mm)
<b>End of measurement 18:05 (24/07/2024)</b>					

## 2 – Activation (“002\_Cu-6Sm\_Activation\_ramp”)

**Temperature:** 28°C – 400°C (T<sub>sp</sub>= 508.9°C) with a ramp of 10°C/minute

**Flow:** 4 mL/min H<sub>2</sub>, 5 mL/min Ar (measured: 4.60 mL/min H<sub>2</sub>, 4.75 mL/min Ar)

8 minutes delay between injection and MS signal

**Pressure:** atmospheric pressure

Pressure “backpressure\_contr2” is set at 0.3 bars (before it was 0.6 bars) so the gases in the massflow are constant

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
2-3	alignment		RT	1		
4	Ref	Cu K-edge	RT	1		Filter3 = 22
5	NEXAFS	Cu K-edge	29°C	1	0.76	Mode 5, <b>start H2 flow</b> , beamstop IN
6	NEXAFS	Cu K-edge	29°C	1	0.77	Mode 5
7	NEXAFS	Cu K-edge	29°C	1	0.77	Mode 5
8	NEXAFS	Cu K-edge	29°C	1	0.77	Mode 5
9	NEXAFS	Cu K-edge	29°C	1	0.77	Mode 5
10	NEXAFS	Cu K-edge	29°C	1	0.77	Mode 5, <b>start</b>

							heating 10°C/min
11	NEXAFS	Cu K-edge	53°C	1	0.78	Mode 5	
12	NEXAFS	Cu K-edge	74°C	1	0.78	Mode 5	
13	NEXAFS	Cu K-edge	105°C	1	0.78	Mode 5	
14	NEXAFS	Cu K-edge	131°C	1	0.78	Mode 5	
15	NEXAFS	Cu K-edge	170°C	1	0.78	Mode 5	
16	NEXAFS	Cu K-edge	193°C	1	0.79	Mode 5	
17	NEXAFS	Cu K-edge	222°C	1	0.79	Mode 5	
18	NEXAFS	Cu K-edge	250°C	1	0.79	Mode 5, starts reduction!	
19	NEXAFS	Cu K-edge	281°C	1	0.79	Mode 5	
20	NEXAFS	Cu K-edge	314°C	1	19:04-19:07	Mode 5 (ends at 327.5°C, starts dwell)	

### 3 – Activation stable (“003\_Cu-6Sm\_Activation\_dwell”)

**Temperature:** 260°C (T<sub>sp</sub>= 327.5°C) starts dwell at 19:07

**Flow:** 4 mL/min H<sub>2</sub>, 5 mL/min Ar (measured: 4.60 mL/min H<sub>2</sub>, 4.75 mL/min Ar)

8 minutes delay between injection and MS signal

**Pressure:** atmospheric pressure

Pressure “backpressure\_contr2” is set at 0.3 bars (before it was 0.6 bars) so the gases in the massflow are constant

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	eiger	18 keV	327.5°C	1	19;13	60 sec
3	Ref	Cu K-edge	327.5°C	1	19;22	Filter3 = 22
4-7	EXAFS	Cu K-edge	327.5°C	1	19;25	Mode 4
8	Ref	Fe K-edge	327.5°C	1	19;57	Filter3 = 19
9-14	NEXAFS	Ce L <sub>3</sub> -edge	327.5°C	1	20;00	Mode 2
15	Ref	Cu K-edge	327.5°C	1	21;12	Filter3 = 22
16						WRONG
17-20	EXAFS	Cu K-edge	327.5°C	1	21;15	Mode 4
21	eiger	18 keV	327.5°C	1	21;50	60 sec

### 4 – Activation cooldown (“004\_Cu-6Sm\_Activation\_CD”)

**Temperature:** From 260 to 200°C (T<sub>sp</sub>= 249.7°C) starts at 22:14. Ends at 22:22.

**Flow:** 9 mL/min Ar (measured: X mL/min Ar)

10°C/min de rampa de enfriamiento

**Pressure:** atmospheric pressure

**Protocol:** eigerloopscan 60 10 (60 scans a 10 segundos por scan)

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	Ref	Cu K-edge	327.5°C	1		Filter3 = 22, Wrong
					22;08	Ar flow to 9 ml/min
2-43	eigerloops	18 keV	327.5-249.7°C	1	22;13	10°C/min

## 5 – Reaction T1=200°C (“005\_Cu-6Sm\_Reaction\_T1”)

**NOTA:** problema con las presiones. El MFC del CO<sub>2</sub> no aguanta más de 4.06 bar de presión. Hasta que se arregle, se va a probar la reacción de metanación (P atmosférica-200-400°C de 50 en 50°C)

**Temperature:** 200°C (T<sub>sp</sub>= 249.7°C) empieza a las 00:25

**Flow:** 4.69 mL/min H<sub>2</sub>\*, 5 mL/min Ar, 1 mL CO<sub>2</sub> (measured: 5.4 mL/min H<sub>2</sub>, 4.75 mL/min Ar, 1.35 mL/min CO<sub>2</sub>)

\*Extrapolando

**Pressure:** 30 bar. Subida de P (10 bar/min) en 50 mL/min Ar (22:26-22:54). Después, cambio a flujos de reacción. **FINALMENTE** Patmosférica.

Pressure “backpressure\_contr2” is set at 0.3 bars (before it was 0.6 bars) so the gases in the massflow are constant

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	eiger	18 keV	249.7°C	1	00;30	60 sec
2	Ref	Cu K-edge	249.7°C	1	00;33	Filter3 = 22
3-6	EXAFS	Cu K-edge	249.7°C	1	00;41	Mode 4
7	Ref	Fe K-edge	249.7°C	1	1;16	Filter3 = 19
8-13	NEXAFS	Ce L <sub>3</sub> -edge	249.7°C	1	1;31	Mode 2
14	Ref	Cu K-edge	249.7°C	1	2;15	Filter3 = 22
15-18	EXAFS	Cu K-edge	249.7°C	1	2;24	Mode 4
19	eiger	18 keV	249.7°C	1	3;02	60 sec

## 6 – Heating from 200 to 250°C (“006\_Cu-6Sm\_Ramp\_T1T2”)

**Temperature:** From 200 to 250°C (T<sub>sp</sub>= 249.7°C to 312.1°C) starts at 3:16. Ends at 3:22.

**Flow:** 4.69 mL/min H<sub>2</sub>\*, 5 mL/min Ar, 1 mL CO<sub>2</sub> (measured: 5.4 mL/min H<sub>2</sub>, 4.75 mL/min Ar, 1.35 mL/min CO<sub>2</sub>)

\*Extrapolando

10°C/min de rampa de calentamiento

**Pressure:** atmospheric pressure

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	Ref	Cu K-edge	249.7°C	1	3;05	Filter3 = 22
2-4	NEXAFS	Cu K-edge	249.7 to 312.1	1	3;16	10°C/min; inicio heating
5-6						WRONG

## 7 – Reaction T<sub>2</sub>=250°C (“007\_Cu-6Sm\_Reaction\_T2”)

**Temperature:** 250°C (T<sub>sp</sub>= 312.1°C) empieza a las 3:22

**Flow:** 4.69 mL/min H<sub>2</sub>\*, 5 mL/min Ar, 1 mL CO<sub>2</sub> (measured: 5.4 mL/min H<sub>2</sub>, 4.75 mL/min Ar, 1.35 mL/min CO<sub>2</sub>)

\*Extrapolando

**Pressure:** P atmosférica

Pressure “backpressure\_contr2” is set at 0.3 bars (before it was 0.6 bars) so the gases in the massflow are constant

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	eiger	18 keV	312.1°C	1	3;33	60 sec
2	Ref	Cu K-edge	312.1°C	1	3;36	Filter3 = 22
3-6	EXAFS	Cu K-edge	312.1°C	1	3;48	Mode 4
7	Ref	Fe K-edge	312.1°C	1	4;21	Filter3 = 19
8-13	NEXAFS	Ce L <sub>3</sub> -edge	312.1°C	1	4;28	Mode 2
14	Ref	Cu K-edge	312.1°C	1	5;15	Filter3 = 22
15-18	EXAFS	Cu K-edge	312.1°C	1	5;23	Mode 4. Beam caído (scan 18, 5;47)
19	eiger	18 keV	312.1°C	1	6;08	60 sec

## 8 – Heating from 250 to 300°C (“008\_Cu-6Sm\_Ramp\_T2T3”)

**Temperature:** From 250 to 300°C (T<sub>sp</sub>= 312.1°C to 379.3°C) starts at 6:24. Ends at 6:30.

**Flow:** 4.69 mL/min H<sub>2</sub>\*, 5 mL/min Ar, 1 mL CO<sub>2</sub> (measured: 5.4 mL/min H<sub>2</sub>, 4.75 mL/min Ar, 1.35 mL/min CO<sub>2</sub>)

\*Extrapolando

10°C/min de rampa de calentamiento

**Pressure:** atmospheric pressure

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	Ref	Cu K-edge	312.1°C	1	6;15	Filter3 = 22
2-4	NEXAFS	Cu K-edge	312.1 to 379.3	1	6;24	10°C/min; inicio heating

## 9 – Reaction T3=300°C (“009\_Cu-6Sm\_Reaction\_T3”)

**Temperature:** 300°C (T<sub>sp</sub>= 379.3°C) empieza a las 6;31

**Flow:** 4.69 mL/min H<sub>2</sub>\*, 5 mL/min Ar, 1 mL CO<sub>2</sub> (measured: 5.4 mL/min H<sub>2</sub>, 4.75 mL/min Ar, 1.35 mL/min CO<sub>2</sub>)

\*Extrapolando

**Pressure:** P atmosférica

Pressure “backpressure\_contr2” is set at 0.3 bars (before it was 0.6 bars) so the gases in the massflow are constant

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	eiger	18 keV	379.3°C	1	6;40	60 sec
2	Ref	Cu K-edge	379.3°C	1	6;43	Filter3 = 22 ??
3	Ref	Cu K-edge	379.3°C	1	6;55	Filter3 = 22 ??
4-5	EXAFS	Cu K-edge	379.3°C	1	7;03	Mode 4 ??
6	Ref	Cu K-edge	379.3°C	1	7;08	Filter3 = 22
7-10	EXAFS	Cu K-edge	379.3°C	1	7;12	Mode 4
11	Ref	Fe K-edge	379.3°C	1	7;41	Filter3 = 19
12-15	NEXAFS	Ce L <sub>3</sub> -edge	379.3°C	1	7;50	Mode 2
16-17	NEXAFS	Ce L <sub>3</sub> -edge	379.3°C	1	8; 32	Mode 2
18	Ref	Cu K-edge	379.3°C	1	8;46	Filter3 = 22
19-22	EXAFS	Cu K-edge	379.3°C	1	8;54	Mode 4
23	eiger	18 keV	379.3°C	1	9;30	60 sec

## 10 – Heating from 300 to 350°C (“010\_Cu-6Sm\_Ramp\_T3T4”)

**Temperature:** From 300 to 350°C (T<sub>sp</sub>= 379.3°C to 444.4°C) starts at 9:4X. Ends at 9:4X.

**Flow:** 4.69 mL/min H<sub>2</sub>\*, 5 mL/min Ar, 1 mL CO<sub>2</sub> (measured: 5.4 mL/min H<sub>2</sub>, 4.75 mL/min Ar, 1.35 mL/min CO<sub>2</sub>)

\*Extrapolando

10°C/min de rampa de calentamiento

**Pressure:** atmospheric pressure

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	Ref	Cu K-edge	379.3°C	1	9:39	Filter3 = 22
2-4	NEXAFS	Cu K-edge	379.3 to 444.4	1	9:42	10°C/min ramp, mode 5
5-6	NEXAFS	Cu K-edge	444.4°C	1	9:54	<b>Change H<sub>2</sub> source</b> (small → big bottle) 16,0 h in MS

## 11 – Reaction T4=350°C (“011\_Cu-6Sm\_Reaction\_T4”)

**Temperature:** 350C (T<sub>sp</sub>= 444.4°C) empieza a las

**Flow:** 4.69 mL/min H<sub>2</sub>\*, 5 mL/min Ar, 1 mL CO<sub>2</sub> (measured: 5.4 mL/min H<sub>2</sub>, 4.75 mL/min Ar, 1.35 mL/min CO<sub>2</sub>)

\*Extrapolando

**Pressure:** P atmosférica

Pressure “backpressure\_contr2” is set at 0.3 bars (before it was 0.6 bars) so the gases in the massflow are constant

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	eiger	18 keV	444.4°C	1	10:03	60 sec
2	Ref	Cu K-edge	444.4°C	1	10:10	Filter3 = 22, <b>BAD!</b>
3-5	EXAFS	Cu K-edge	444.4°C	1	10:14	Mode 4, ABORT
6-9	EXAFS	Cu K-edge	444.4°C	1	10:26	Mode 4, beamdump in last scan (point 122)
10	Ref	Cu K-edge	444.4°C	1	11:35	After injection, <b>BAD!</b>
11-14	EXAFS	Cu K-edge	444.4°C	1	11:38	
15	Ref	Fe K-edge	444.4°C	1	12:08	Filter3 = 19, <b>BAD!</b>
16-21	NEXAFS	Ce L <sub>3</sub> -edge	444.4°C	1	12:14	Mode 2, ABORT 21
22-23	NEXAFS	Ce L <sub>3</sub> -edge	444.4°C	1	13:01	Mode 2
24	Ref	Cu K-edge	444.4°C	1	13:15	Filter3 = 22, <b>BAD!</b>
27	Ref	Cu K-edge	444.4°C	1	13:40	Filter3 = 22
28-31	EXAFS	Cu K-edge	444.4°C	1	13:44	Mode 4
32	eiger	18 keV	444.4°C	1	14:23	60 sec

- BAD** means that the ionization chamber (“ionch1”) was frozen and hat to be restarted with “Keithley\_reset\_6485”



## 12 – Heating from 350 to 400°C (“012\_Cu-6Sm\_Ramp\_T4T5”)

**Temperature:** From 350 to 400°C (T<sub>sp</sub>= 444.4°C to 508.9°C) starts at 14:XX. Ends at 15:XX.

**Flow:** 4.69 mL/min H<sub>2</sub>\*, 5 mL/min Ar, 1 mL CO<sub>2</sub> (measured: 5.4 mL/min H<sub>2</sub>, 4.75 mL/min Ar, 1.35 mL/min CO<sub>2</sub>)

\*Extrapolando

10°C/min de rampa de calentamiento

**Pressure:** atmospheric pressure

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	Ref	Cu K-edge	444.4°C	1	14:45	Filter3 = 22
2-5	NEXAFS	Cu K-edge	444.4 to 508.9	1	14:57	10°C/min ramp, mode 5
				1		

## 13 – Reaction T4=400°C (“013\_Cu-6Sm\_Reaction\_T5”)

**Temperature:** 400C (T<sub>sp</sub>= 508.9°C) empieza a las 15:08h

**Flow:** 4.69 mL/min H<sub>2</sub>\*, 5 mL/min Ar, 1 mL CO<sub>2</sub> (measured: 5.4 mL/min H<sub>2</sub>, 4.75 mL/min Ar, 1.35 mL/min CO<sub>2</sub>)

\*Extrapolando

**Pressure:** P atmosférica

Pressure “backpressure\_contr2” is set at 0.3 bars (before it was 0.6 bars) so the gases in the massflow are constant

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	eiger	18 keV	508.9°C	1	15:13	60 sec
2	Ref	Cu K-edge	508.9°C	1	15:19	Filter3 = 22,
3-6	EXAFS	Cu K-edge	508.9°C	1	15:24	Mode 4
7	Ref	Fe K-edge	508.9°C	1	15:54	Filter3 = 19
8-13	NEXAFS	Ce L <sub>3</sub> -edge	508.9°C	1	16:02	Mode 2
14-15	NEXAFS	Ce L <sub>3</sub> -edge				Extra, don't use
16	Ref	Cu K-edge	508.9°C	1	16:50	Filter3 = 22
17-20	EXAFS	Cu K-edge	508.9°C	1	16:53	Mode 4
21	eiger	18 keV	508.9°C	1	17:29	60 sec

## 14 – Cool down from 400 to 20°C (“014\_Cu-6Sm\_CD\_T5RT”)

**Temperature:** From 400 to 20°C (T<sub>sp</sub>= 25°C to 25°C) starts at 17:37. Ends at 19:29h.

**Flow:** 5 mL/min Ar (measured: 4.75 mL/min Ar)

10°C/min de rampa de calentamiento

**Pressure:** atmospheric pressure

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
S1: 0-260	eiger	18keV	32	1	17;36	10 sec each scan; 17;37 starts cooldown

## 15 – Post Mortem at RT (“015\_Cu-6Sm\_PM\_RT”)

**Temperature:** RT starts at 17:37. Ends at 00:28.

**Flow:** 5 mL/min Ar (measured: 4.75 mL/min Ar) IMPORTANT in scan 3 there is no more Ar flow. Dirk was working with the gases so we have our reactor in Ar but no flow.

**Pressure:** atmospheric pressure

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	eiger	18 keV	30	1	18;34	60 sec
2	Ref	Mo k edge	30	1	18;37	Filter3 = 14
3-8	NEXAFS	Zr k edge	30	1	18;42	3 XANES, Mode 2
9	Ref	Cu	30	1	19;21	Filter3 = 22
9	Ref	Cu	30	1		Filter3 = 22
10-12	EXAFS	Cu K-edge				Measurement Stopped by mistake
13-24	EXAFS	Cu K-edge	30	1	19;31	3 EXAFS, Mode 4
25	Ref	Fe K-edge	30	1	21;00	Filter3 = 19
26-27	NEXAFS	Ce L <sub>3</sub> -edge	30	1	21;06	3 NEXAFS, Mode 2
28-29						Checking ropi and roll counts are low for Ce
30-31	NEXAFS	Ce L <sub>3</sub> -edge	30	1	21;26	3 NEXAFS, Mode 2, ERROR Lox counts again
32-37	NEXAFS	Ce L <sub>3</sub> -edge, Fe ref				Several tests
45-50	NEXAFS	Ce L <sub>3</sub> -edge	30	1	23;11	3 NEXAFS, Mode 2
51	Ref	Fe K-edge	30	1	23;55	Filter3 = 19
52-55	NEXAFS	Sm L <sub>3</sub> -edge	30	1	00;01	3 NEXAFS, Mode 2

2024\_07\_25

## Sample: #2 CuFe-6Sm

**Loading:** 2-3 mg of sample in the capillary

**Sample:** 60%wt CuO, 40%wt 6Sm (calcinated at 500°C, 4%at. Zr, 6%at Sm, 90%at Ce)

### 1 – Pristine RT (“001\_CuFe-6Sm\_Pristine\_RT”)

**Temperature:** 28°C

**Flow:** set-point at 5 mL/min Ar (after 3am, 0 mL/min)

**roi:** 1: Ce; 2: Sm; 3: Cu; 4:Zr, 5:Fe

mz = - 1.68

**Pressure:** atmospheric pressure

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	eiger	18 keV	30	1	00;55	60 sec
2	Ref	Cu	30	1	1;03	Filter3 = 22
3-5	EXAFS	Cu K-edge	30	1	1;06	Mode 4. Alinear
6	ropi		30	1	1;11	Para Cu
7	roll		30	1	1;15	Para Cu
8	Ref	Fe	30	1	1;17	Filter3 = 19, <b>BAD</b>
9	ropi		30	1	1;22	For Fe
10	roll		30	1	1;27	For Fe
11	Ref	Fe K-edge	30	1	1;29	Filter = 19
12-15	EXAFS	Fe k-edge	30	1	1;33	Mode 4, roi 5
16	Ref	Cu	30	1	2;08	Filter3 = 22
17-20	EXAFS	Cu K-edge	30	1	2;13	Mode 4
21	Ref	Fe K-edge	30	1	2;42	Filter3 = 19
22-27	NEXAFS	Ce L <sub>3</sub> -edge	30	1	2;48	Mode 2
28	Ref	Fe K-edge	30	1	3;41	Filter3 = 19

### 2 – Activation (“002\_CuFe-6Sm\_Activation\_ramp”)

**Temperature:** 28°C – 260°C (T<sub>sp</sub>= 327.5°C) with a ramp of 10°C/minute. Starts at 3:52. Ends at 4:21

**Flow:** 4 mL/min H<sub>2</sub>, 5 mL/min Ar (measured: 4.60 mL/min H<sub>2</sub>, 4.75 mL/min Ar)

8 minutes delay between injection and MS signal

**Pressure:** atmospheric pressure

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1-10	NEXAFS	Fe K-edge	28 to 327.5°C	1	3;52	10°C/min ramp, mode 5

### 3 – Activation stable (“003\_CuFe-6Sm\_Activation\_dwell”)

**Temperature:** 260°C (T<sub>sp</sub>= 327.5°C) starts dwell at 4:22

**Flow:** 4 mL/min H<sub>2</sub>, 5 mL/min Ar (measured: 4.60 mL/min H<sub>2</sub>, 4.75 mL/min Ar)

8 minutes delay between injection and MS signal

**Pressure:** atmospheric pressure

#	Type	Scan info	Temp(T <sub>sp</sub> )	Pres (bar)	Time	Comments
1	Ref	Cu	327.5	1	4;25	Filter3 = 22
2-5	EXAFS	Cu K-edge	327.5	1	4;28	Mode 4
6	Ref	Fe	327.5	1	4;58	Filter3 = 19
7-14	NEXAFS	Fe k-edge	327.5	1	5;05	Mode 2. scan18 (beam off 5;45-6;05)
15	Ref	Fe	327.5	1	6;09	Filter3 = 19. <b>BAD</b>
16	ropi		327.5	1	6;12	Para Cu. <b>BAD</b>
17	Ref	Fe	327.5	1	6;17	Filter3 = 19
18-19	NEXAFS	Fe k-edge	327.5	1	6;20	Mode 2. <b>BAD</b>
20-21	NEXAFS	Ce L <sub>3</sub> -edge	327.5	1	6;35	Mode 2
22	ropi		327.5	1	6;45	Para Ce
23	roll		327.5	1	6;47	Para Ce
24-25	NEXAFS	Ce L <sub>3</sub> -edge	327.5	1	6;49	Mode 2. <b>Pocas cuentas</b>
26	eiger	18 keV	327.5	1		60 sec
27	eiger	18 keV	327.5	1	7;16	60 sec; check evol.
28	eiger	18 keV	327.5	1	7;24	60 sec; check evol.

### 4 – Activation cooldown (“004\_CuFe-6Sm\_Activation\_CD”)

**Protocol:** eigerloopscan 60 10 (60 scans a 10 segundos por scan)